

MEWS Development Program

TEST CASE DOCUMENT — Critical Alert Triggering and Display

Test Case ID	DR-TC-002	Test Type	Functional / Alerting
Feature	MEWS Alerting	Execution Status	Not Executed
Persona	Clinical Staff	Priority	Critical
Test Basis	Dry Run Test Case	Environment	To be completed during execution

1. Test Scenario

Critical Alert Triggering and Display

A patient is admitted with vital signs simulating a MEWS score of ≥ 5 . The test verifies that the system calculates the critical score and triggers the corresponding critical/red alert for clinical staff.

2. Preconditions

- Patient record is available in the system.
- Patient is associated with a connected monitoring device.
- Clinical staff user is logged in and monitoring the dashboard.
- Test data represents a critical condition with MEWS ≥ 5 .

3. Test Data

Patient / Test Subject	Patient A (anonymized)
Heart Rate	125 bpm
Respiratory Rate	32 breaths/min
Blood Pressure	130/85 mmHg
MEWS Condition	MEWS ≥ 5
Expected Alert	Critical / Red

4. Test Steps

1	Monitor the dashboard when the test patient is admitted / becomes active.
2	Confirm the patient vital-sign data represents the critical test condition.
3	Observe the MEWS calculation and patient status displayed by the system.
4	Observe the alert mechanism presented to the clinical staff user.
5	Record the actual system behavior and attach supporting evidence, if available.

5. Expected Result

The system automatically calculates **MEWS ≥ 5** , triggers a **Critical / Red Alert** to the clinical staff user, and prominently displays the critical status on the dashboard.

6. Actual Result

Actual Result	
---------------	--

7. Evidence / Screenshots

Evidence	Attach screenshot, test log, or execution record here.
-----------------	--

8. Result & Review

Status	☐ Pass ☐ Fail ☐ Blocked ☐ Not Executed
Defect / Issue ID	_____
Tester	_____
Execution Date	_____
Reviewer	_____

Portfolio sample. This document is based on the supplied Dry Run Test Case structure and source scenario. It is anonymized and does not represent a controlled QMS record or claim an executed test result.